

Stroke Volume Variation

“Can We Use Fluid to Improve Hemodynamics?”

Introduction:

In the quest to achieve optimal oxygen delivery (DO_2), clinicians are often forced to use imprecise, non-specific information to guide their therapy. Traditional hemodynamic monitoring parameters (HR, MAP, CVP, and PAOP) are often insensitive and sometimes misleading in the assessment of circulating blood volume. However, the appropriateness of their interventions is often crucial to avoid the deleterious effects of over-, under-, or inappropriate resuscitation. Volume is one of the first therapeutic interventions selected when optimizing DO_2 . Often times the choice to intervene using fluid is accompanied by the difficult questions, “Can using fluid improve hemodynamics?” and, “Is it the appropriate intervention?” Stroke volume variation (SVV) as available on the FloTrac system may help answer these questions.


What Causes Stroke Volume Variation?

Stroke volume variation is a naturally occurring phenomenon in which the arterial pulse pressure falls during inspiration and rises during expiration due to changes in intra-thoracic pressure secondary to negative pressure ventilation (spontaneously breathing). Variations over 10mmHg have been referred to as pulsus paradoxus. The normal range of variation in spontaneously breathing patients has been reported between 5-10mmHg.

Reverse pulsus paradoxus is the same phenomenon with controlled mechanical ventilation, however, in reverse. Arterial pressure rises during inspiration and falls during expiration due to changes in intra-thoracic pressure secondary to positive pressure ventilation. In addition to reverse pulsus paradoxus, it has also been referred to as paradoxical pulsus, respiratory paradox, systolic pressure variation and pulse pressure variation. Traditionally SVV is calculated by taking the $SV_{max} - SV_{min} / SV_{mean}$ over a respiratory cycle or other period of time.

SVV and Assessing Fluid Response

SVV and its comparable measurement, pulse pressure variation (PPV), are not indicators of actual preload but of relative preload responsiveness. SVV has been shown to have a very high sensitivity and specificity when compared to traditional indicators of volume status (HR, MAP, CVP, PAD, PAOP), and their ability to determine fluid responsiveness. The following table of studies demonstrates SVV sensitivity and specificity in predicting fluid responsiveness against a specified infused volume and defined criteria for a fluid responder.



Study	Patients	Volume	Tidal Volume ml/Kg	Parameters Tested (Artery)	R ²	Def. of Responder	Sensitivity	Specificity
Michard ²	Sepsis	500 ml	8 to 12	ΔPP (R or F)	0.85	ΔCO ₂ ≥ 15%	94	96
Berkenstadt, et al. ¹	Neuro Surgery	100 ml	10	ΔSVV	0.53	ΔSV ≥ 5%	79	93
Reuter, et al. ³	Cardiac	10 x BMI	10	ΔSVV	0.64	ΔSV ≥ 5%	79	85

How Can I Use SVV?

Normal SVV values are less than 10-15% on controlled mechanical ventilation. The figures to the right demonstrate using SVV as a guide for volume resuscitation with a goal SVV of <13%. SVV increased to 19% with a stroke volume (SV) of 43 ml/beat, blood and saline were given to obtain a SVV of 6% and a SV of 58 ml/beat.

What are the Limitations and Effects of Other Therapies on SVV?

• Mechanical Ventilation

Currently, literature supports the use of SVV only on patients who are 100% mechanically (control mode) ventilated with tidal volumes of more than 8cc/kg and fixed respiratory rates.

• Spontaneous Ventilation

Currently, literature does not support the use of SVV with patients who are spontaneously breathing due to the irregular nature of rate and tidal volumes.

• Arrhythmias

Arrhythmias can dramatically affect SVV values. Thus, the utility of SVV as a guide for volume resuscitation is greatest in absence of arrhythmias.

• PEEP

Increasing levels of positive end expiratory pressure (PEEP) may cause an increase in SVV, the effects of which may be corrected by additional volume resuscitation if warranted.

• Vascular Tone

The effects of vasodilatation therapy may increase SVV and should be considered before treatment with additional volume.

Summary



When used within its limitations SVV is a sensitive tool that can be used to guide the appropriate management of the patient's preload to achieve optimal DO₂ and answer the question "Can we use fluid to improve hemodynamics?" SVV is an available parameter with the FloTrac sensor and the Vigileo monitor.

*Continuous SVR requires CVP to be supplied to Vigileo monitor from a bedside monitor.

Prepared by: J. Frazier RN, RRT, Clinical Marketing, Edwards Lifesciences

References

- Berkenstadt H, et al. Stroke volume variation as a predictor of fluid responsiveness in patients undergoing brain surgery. *Anesth Analg* 2001;92:984-989.
- Michard F. Changes in arterial pressure during mechanical ventilation. *Anesthesiology* 2005; 103:419-428.
- Reuter DA, et al. Usefulness of left ventricular stroke volume variation to assess fluid responsiveness in patients with reduced cardiac function. *Crit. Care Med* 2003; 31:1300-404.



Rx only. See instructions for use for full prescribing information.

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

Edwards and Vigileo are trademarks of Edwards Lifesciences Corporation. Edwards Lifesciences and the stylized E logo and FloTrac, are trademarks of Edwards Lifesciences Corporation and are registered in the United States Patent and Trademark Office.

© 2007 Edwards Lifesciences LLC
All rights reserved. AR02897



Edwards Lifesciences